

Patient _____

Doctor _____ Date _____

VARISAN[®] TOP

CUSTOM MADE - CIRCULAR KNIT



CIRCULAR KNIT



CUSTOM MADE



- | | | | | | |
|--|---------------------------------|-------------------------------------|--------------------------------|--------------------------------|--------------------------------|
| ☐ Ccl. 1 18 - 21 mmHg | ☐ COTTON | ☐ OPEN TOE | ☐ PIECE | ☐ Beige | ☐ Black |
| ☐ Ccl. 2 23 - 32 mmHg | | ☐ CLOSED TOE | ☐ PAIR | | |
| ☐ Ccl. 3 34 - 46 mmHg | | | | | |


POINT FOR MEASUREMENTS:

Patient lay on the couch (circumferences)

- Point cA: circumference at the base of the toes
 Point cY: circumference from the lower margin of the heel to the instep
 Point cB: circumference of the narrowest point found above the malleoli
 Point cC: circumference at the widest point of the calf
 Point cD: circumference at the head of the fibula
 Point cF: midway circumference between point E and point G

Patient standing (circumferences and lengths)

- Point cG: circumference at the widest point of the thigh, 2-3 cm under the groin
 Point K: identify the lengths: front (from the groin to the waist), back (from the gluteal fold to the waist)
 Point cH: circumference at the widest point of the pelvis
 Point cT: circumference at the waistline



☐ KNEE-HIGH STOCKING	
☐ SELF-SUPPORTING KNEE-HIGH STOCKING	
☐ KNEE-HIGH STOCKING WITH ZIP	
☐ SIDE ZIP	☐ RH ☐ LH
☐ BACK ZIP	
☐ FULL-LENGTH STOCKING	
☐ SELF-SUPPORTING STOCKING	
☐ SINGLE-LEG TIGHTS	☐ RH ☐ LH
☐ TIGHTS	☐ MAN ☐ WOMAN
☐ MATERNITY TIGHTS	

VARITEC[®] STOCKING AID

☐ Ref. 9705008

VARISLIDE[®] SLIDER

☐ Short ready-up
Ref. 970510C

☐ Long ready-up
Ref. 970510L